




Optimice®

- Highest housing density in the industry
- Fully protects research animals
- No motors, blowers, noise, or vibration
- Prevents release of odors and allergens in the animal room
- Creates a variable-free environment
- Lightweight, easy and low cost to maintain
- Automatic watering options
- Single-use caging available; interchangeable in Optimice racks
- Many customized accessories available

Specifications

Rack Capacity: 100 cages
Cage Floor Space: 75 in² (484 cm²)
Cage Interior Height: 5 inches (13 cm)
Water Capacity: 13 oz (400 mL) bottle or external bag
Feeder Capacity: 17 oz (510 mL)
Rack Casters: 4 swivel, 2 locking

	Optimice Rack					
	Dimensions <i>Inches (mm)</i>				Weight <i>Pounds (kg)</i>	
	Depth	Width	Height	Auto-Water Height	Basic	Auto-Water
<i>Full rack with card holders</i>	40.5 (1029)	40.5 (1029)	76.5 (1943)	78 (1981)	613 (278)	659 (299)
<i>Full rack no card holders</i>	39.7 (1001)	39.7 (1001)			604 (274)	650 (295)
<i>Empty rack</i>	33.7 (856)	39.7 (1001)			326 (148)	366 (166)

All Animal Care Systems products and accessories are listed and described at animalcaresystems.com.