



SAY HELLO
TO EXPERTISE

MD10
(Magnetic Stirrer)

MD10

(Magnetic Stirrer)

When looking for magnetic stirrers, users look for a device which delivers reliable performance and years of durability, but most manufacturers tend to oversell and under deliver which results in a frustrating experience for the user. Magnetic stirrers come with promises of exceptional performance, quality components and high durability, but in just a short time, the reality becomes apparent. Users face grave issues like the top getting corroded or the motor is not powerful enough to stir large quantities of solutions, its speed remains low etc.

With detailed insights about the requirements and needs of users and their feedback of existing stirrers, the R&D department at Stirrlab came up with the design of iStir MD10 magnetic stirrer to develop a product that delivers high performance without compromising on flexibility or durability. With its stainless steel top surface, the device is extremely resistant to corrosion while its long lasting DC motor runs at top speed without needing much of a maintenance. The same motor is capable enough to stir large quantities of solutions at once for higher productivity and less wastage of time and effort.



Product Features



Stirring Capacity
Up To 10 Litres



Variable Speed Setting
From 200-2200 RPM With
Increments Of 10 Rpm




Variable Timer Setting
From 1 To 999 Minutes
Along With An Infinite Mode.



Powerful DC Motor
For Quite & Long Life
Performance

Microprocessor Controlled
For Efficient Performance
Even Under High Loads.



Digital Display Allows The
User To Gather All The Relevant
Information In A Single Glance.



Pulse Mode Function With
Run-pause-run Function For
Thorough Mixing Of Solutions.



No corrosive SS top
increases life of
instrument



Product Features



Closed Design Allows
Easy Cleaning And
Protects The Equipment
From Accidental Spills.



Non-slip Strength Rubber
Feet For Maximum Stability
Even At High-speeds.



Compact Footprint Saves
Valuable Bench Space And
Also Makes It Suitable For
Use In Hoods And Cold Room.

IP42 Class Compliant
Rating For Unmatched
Durability.



Safe Low
Voltage Device



Specifications

Specifications	MD10
Capacity	10 Litres
Variable Speed	200-2200 RPM
Run Time	1-999 mins & infinite mode
Pulse Mode	Yes, (Run - pause - Run Mode)
Digital Display	Yes
Weight	1.4 kg
Dimension	156 x 248 x 104 mm
Protection class	IP 42
Setup plate Material	Plastic Integrated SS Top Plate

Ordering Information

MD10		
Product Name	Cat #	Details
MD10	HS20-2314 EU/US/AU/UK	Magnetic Stirrer

Standard Accessories	
Part Name	Cat #
Stir Bar	B1



www.labquipindoprima.co.id